



Laerdal
helping save lives



Nursing Anne

Effective Clinical Skills Training

The Nursing Anne manikin is educationally effective for clinical training, targeting key skills from basic patient care to advanced nursing competencies. Nursing Anne paired with the SimPad PLUS System provides an economical and relevant training platform for women's health, wound assessment and general patient care.

Nursing Anne

Patient Care:

- Head with anatomical landmarks, trachea, and esophagus, along with simulated lungs and stomach, allow the practice of many procedures
- Irrigation of the eye and ear (simulated)
- Application/installation of medications in the eye, ear and nose including nasal packing
- Mouth and denture care procedures
- Insertion and suctioning of oropharyngeal and nasopharyngeal airways
- Insertion, securing, and care of endotracheal tubes
- Tracheostomy care and tracheal suctioning
- Various oxygen delivery procedures
- NG tube insertion, care, medication administration, and removal (Updated 2019)
- Gastric lavage and gavage
- Full range of motion for realistic patient handling
- Fingers and toes are spread to allow bandaging
- Interchangeable stomas depict colostomy, ileostomy and suprapubic cystostomy
- Colostomy may be irrigated and will retain an indwelling catheter

Interchangeable Male and Female Genitalia:

- Complete urinary catheterization
- Can be attached to urinary and colon reservoirs via connector valves
- Female genitalia capable of vaginal douching
- Will retain indwelling or straight catheter
- Enema procedures may be performed using fluid for realistic return
- When used with manikin, fluid may be used for realistic return
- Urinary valves give the natural resistance felt when catheterizing
- Anal valves simulate the internal anal sphincter

Circulatory Skills and IV Drug Administration:

- Articulating IV training arm with replaceable skin and infusible vein system allows peripheral intravenous therapy and site care
- Venipuncture possible in the antecubital fossa and dorsum of the hand
- Accessible veins include median, basilic and cephalic
- Subclavian, jejunostomy and Hickman catheter openings (placement only)
- Deltoid, dorsogluteal, and vastus lateralis IM injections possible

Blood Pressure / Pulses:

- Articulating blood pressure arm for auscultated and palpated blood pressure simulation
- Korotkoff sounds synchronized with ECG
- Korotkoff sounds volume control in 10 steps, 0-9
- Systolic and diastolic pressure may be set individually in steps of 2 mmHG
- Systolic 0-300 mmHG, Diastolic 0-200 mmHG
- Auscultative Gap, with on/off feature
- Pressure accuracy +/- 4 mmHG
- Calibrate function to adjust pressure sensor and cuff gauge
- Palpated pulse simulator controls brachial and radial pulses
- Pulses only active when palpated
- Pulses synchronized with programmable ECG
- Pulse strengths dependent on Blood Pressure
- Manually generated carotid pulse

Sounds - SimPad PLUS System:

- Heart sounds synchronized with ECG
- Auscultated lung sounds synchronized with breathing, 0 - 60 BPM
- Individual lung sound selection
- Normal or abnormal bowel sounds
- Vocal sounds: Computer-generated sounds, recorded vocal sounds and real-time voice input via headset
- User generated vocal sounds

SimPad PLUS System Function:

- Time stamped activities, vital signs, and notes are captured in the data log recording alongside the time-stamped data log
- The SimPad PLUS log files can be viewed on the device or on a PC for post-simulation reflection and debriefing
- Upload self-authored scenarios and Themes, or download pre-programmed scenarios
- Handheld, intuitive touchscreen remote for easy 'pick up and play' experience
- Mobile - teach anywhere
- Operate on-the-fly or utilize scenarios and Themes for consistent simulation training

Patient Monitor (Optional):

- Touch screen simulated patient monitor provides concise clinical feedback for physiological parameters
- The monitor's color screen is configurable and provides multiple simulated parameters, each presenting multi-level alarms
- Simulated parameters include HR, ECG, SpO2, BP, RR, Temperature, and etCO2

ORDERING INFORMATION

325-05050 Nursing Anne Light Skin

325-05050T Nursing Anne Tan Skin

325-05050B Nursing Anne Brown Skin

Includes: Adult Manikin, Hospital Gown, Airway Lubricant, Assembly Tool Kit and Directions for Use.

OPERATION AND PATIENT MONITORING OPTIONS

204-30001 SimPad PLUS System

Includes: SimPad PLUS Handheld Remote Control, SimPad PLUS Link Box, Lithium Ion Battery, AC Adapter and Power Cord, USB Cable, Ethernet Cable, Wrist Strap, SimPad PLUS Protective Sleeve, SimPad PLUS Nametags, Manikin Strap, Headset, Microphone and Quick Start Guide.

204-50150 LLEAP for SimPad PLUS Software License

400-10201 Laptop Patient Monitor

400-09201 Tablet Patient Monitor

400-09501 Rugged Patient Monitor

400-29301 All in One Panel PC Patient Monitor

ACCESSORIES AND CONSUMABLES

325-00450, 325-00450T, 325-00450B Fundus Skills and Assessment Module

325-00550, 325-00550T, 325-00550B Wound Care and Assessment Set

325-00650, 325-00650T, 325-00650B Mastectomy Module

205-10150, 205-10150T, 205-10150B Adult IO Leg

325-0075 Breast Exam Module

Operation

Nursing Anne can be operated by a SimPad/SimPad PLUS remote and can utilize a simulated patient monitor. Incorporate SimCapture for a complete debriefing solution.



Optimize your investment

Laerdal offers comprehensive Educational and Technical Services designed to support healthcare professionals successfully implement and sustain simulation training programs. Our team of educators and technicians have years of industry related experience, with the ability to provide a wide range of options to help you maintain your simulation equipment and gain the critical skills necessary to optimize performance and learning outcomes.



Laerdal
helping save lives